

Tutorialal:

Data Analysis for Multi-channel EEG Recordings during A Sustained-attention Driving Task

I. Overview


I.1. Software, toolbox, and environment:

- MATLAB R2007a
- EEGLAB v. 5.03
- Running under workstation (Linux): BRC server 140.113.34.3

I.2 Pre-requirements


- Critical reasoning and thinking
- Being curious and skeptical about the (discovered) results
- Basic MATLAB script coding ability
- EEGLAB usage
- English communication skill
- C language programming (suggestion)

I.3 Flowchart



- Do NOT change the sequence of these procedures unless without previous notification.
 - If the procedure can be done later, it will be specified in this documentation.

II. Integrating behavioral data into EEG signals



- MATLAB scripts
 - `eventnewDN.m`
 - `DR_plot_traj.m`
- Procedures
 - Integrate “log” file (from stimulus computer) and “ev2” (from Neuroscan) file
 - ➔ “event” file (`sxx_yymmdd(MN)_event.txt`)
 - ◆ Use MATLAB script `eventnewDN.m` to integrate “log” file and “ev2” file

```
1 tic
2 subj = 's44';
3 ExpDate = '070205';
4 MainDIR = ['~/liang/' subj '_' ExpDate '/'];
5 % subj = 's31';
6 % ExpDate = '061020';
7 % MainDIR = ['~/liang/' subj '_' ExpDate '/'];
8 SR = 500; %sampling rate of EEG when RECORDING
9
10 A=[];
11 B=[];
12 C=[];
13 D=[];
14 E=[];
15 F=[];
16 A3=[];
17 A4=[];
18 A = importdata([MainDIR subj '_' ExpDate '.ev2']);
19 % A = importdata([MainDIR 's42nm-070105.ev2']);
20 if isstruct(A)
21     A = A.data;
22 end
23 % B = load([MainDIR subj '_' ExpDate '_log.txt']);
24 B = load([MainDIR subj '_' ExpDate '.txt']);
25 % B = load([MainDIR subj '_' ExpDate '.txt']);
26 % B = load([MainDIR 's42nm-070105.txt']);
27
28 A = sortrows(A, 6);
29 for i = size(A, 1) : -1 : 2
30     if A(i, 6) == A(i - 1, 6)
31         A(i, 7) = 1;
32     end
33 end
34 idx = find(A(:, 7) == 1); A(idx, :) = [];
35 A(:, 7) = [];
```

- ◆ Variables: checking before running
 - subj: subject code
 - ExpDate: date of experiment
 - MainDIR: directory of datasets
 - A: contents of ev2 file
 - B: contents of log file
- Fix response offset (by Darkflame 昂穎 or RichYe 人慈)
- Use MATLAB script DR_plot_traj.m to check if some events miss in the event file
 - Note: this code need DR_check_incomplete.m to check if incomplete trials (trials missing some of the events) exist

```

1 % Plot deviation, for Drowsiness
2 % Required input: event file (subj_ExpDate_event.txt)
3 % Output: subj_ExpDate_driving_err.mat, containing moving averaged trajectory and latency
4 % by poem, 20080430
5
6 % clear all;
7 FilePath = '';
8 if isunix
9     mypath;
10    % FilePath = '~/Drowsiness/';
11    % FilePath = '~/liang/';
12    SL = '/';
13 end
14 SET = {
15     's48_080501n';
16     's48_080516m';
17     's22_080513m';
18     's41_080520m';
19     's49_080527n';
20     's22_080529n';
21     's41_080530n';
22     's49_080602m';
23     's50_080725n';
24     's50_080731m';
25     's54_081209n';
26     's22_090120n';
27 };
28 SR = 500; %DOWN-SAMPLED sampling rate. 250 => 250 Hz
29 % window = 90; %90-sec window
30 % stepping = 2; %2-sec stepping
31 COLOR = {
32     [.5 .5 .5]; %raw traj.
33     [0 0 0]; %moving avg. traj.
34     [1 0 0]; %dev_on (left)
35     [1 0 0]; %dev_on (right)
36     [0 1 0]; %act_on
37     [0 0 1]; %act_off (correct)
38     [.5 .5 .5]; %act_off (wrong)
39 };
40 FilePathOld = FilePath;
41
42 for i = 1 : size(SET, 1)
43     FilePath = FilePathOld;
44     Set = SET{i};
45     [subj ExpDate] = strtok(Set, '_'); ExpDate = ExpDate(2 : end);
46     FilePath = [FilePath Set SL];
47
48     %load modified event file
49     try
50         evt = load([FilePath 'single/' Set '_event_new.txt']); %sxx_yymmdd_event_new.txt
51         evt = load([FilePath Set '_event.txt']); %sxx_yymmdd_event_new.txt
52     catch
53         error('Fail to load modified event file');
54     end
55     evt(:, 2) = mod(evt(:, 2), 1000);
56     %fix the first and the last events
57     if evt(1, 2) > 250 %1st point
58         evt(1, 2) = evt(2, 2);
59     end
60     if evt(end, 2) > 250 & (evt(end, 2) ~= 254 & evt(end, 2) ~= 255) %last point
61         evt(end, 2) = evt(end - 1, 2);
62     end
63     %fix index of event
64     if ~isequal(evt(:, 1), (1 : size(evt(:, 1), 1)))
65         evt(:, 1) = (1 : size(evt(:, 1), 1));
66     end
67     %extract event mark
68     evt_tmp = find(evt(:, 2) >= 250); %evt: event file
69     evt_mark = evt(evt_tmp, [1 2 5]);
70     help DR_check_incomplete;
71     [evt_mark evt_tmp] = DR_check_incomplete(evt_mark, 2);
72     if ~isempty(evt_tmp)
73         fprintf('Manually remove incomplete trials: ');
74         keyboard;
75     end
76     %after this step, evt_mark should be [(251/252) - 253 - (254/255)] + [(251/252) - 253 - (254/255)] ...

```

◆ Variables: check before running

- FilePath: directory of datasets
- SET: trajectories of datasets to be plot (multiple datasets allowed)


- evt: name of event file (notice the format)
 - ◆ If the program does not stop in line 74 shown in the above figure (no missing events) → begin pre-processing
- Else,
- Write down the values in the variable `evt_tmp` (see MATLAB workspace). These values are the indices of rows in the variable `evt_mark` with missing events

| | 1 | 2 | 3 | 4 | 5 | 6 |
|-----|-------|-----|--------|---|---|---|
| 527 | 93120 | 253 | 779187 | 2 | 1 | |
| 528 | 93258 | 254 | 780340 | 3 | 1 | |
| 529 | 93580 | 252 | 783031 | 1 | 1 | |
| 530 | 93685 | 253 | 783909 | 2 | 1 | |
| 531 | 93686 | 254 | 783917 | 3 | 1 | |
| 532 | 94369 | 251 | 789627 | 1 | 1 | |
| 533 | 94416 | 253 | 790020 | 2 | 1 | |
| 534 | 94532 | 254 | 790989 | 3 | 1 | |
| 535 | 94874 | 252 | 793847 | 1 | 1 | |
| 536 | 94923 | 253 | 794256 | 2 | 1 | |
| 537 | 94924 | 254 | 794264 | 3 | 1 | |
| 538 | 95514 | 252 | 799195 | 1 | 1 | |

- Run the remaining program
- Check the rows in `evt_mark` around the indices (in `evt_tmp`). The 3rd column is the latencies in the produced event file (the same as the 5th column in the variable `A3`)


| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 |
|-------|-------|-----|---|---|--------|---------|-----|-----|---------|----|
| 93682 | 93682 | 233 | 0 | 0 | 783884 | 0 | 0 | 0 | 1566610 | |
| 93683 | 93683 | 233 | 0 | 0 | 783892 | 0 | 0 | 0 | 1566625 | |
| 93684 | 93684 | 233 | 0 | 0 | 783900 | 0 | NaN | 0 | 1566641 | |
| 93685 | 93685 | 253 | 0 | 0 | 783909 | -0.7143 | NaN | NaN | 1566656 | |
| 93686 | 93686 | 233 | 0 | 0 | 783917 | -0.7143 | 0 | NaN | 1566672 | |
| 93687 | 93687 | 233 | 0 | 0 | 783926 | -0.7143 | 0 | 0 | 1566688 | |
| 93688 | 93688 | 233 | 0 | 0 | 783934 | -1.4286 | 0 | 0 | 1566703 | |
| 93689 | 93689 | 233 | 0 | 0 | 783942 | -1.4286 | 0 | 0 | 1566719 | |
| 93690 | 93690 | 233 | 0 | 0 | 783951 | -1.4286 | -1 | 0 | 1566735 | |
| 93691 | 93691 | 232 | 0 | 0 | 783959 | -1.4286 | 0 | -1 | 1566750 | |
| 93692 | 93692 | 232 | 0 | 0 | 783967 | -2.1429 | 0 | 0 | 1566766 | |
| 93693 | 93693 | 232 | 0 | 0 | 783976 | -2.1429 | -1 | 0 | 1566797 | |
| 93694 | 93694 | 231 | 0 | 0 | 783984 | -2.1429 | 0 | -1 | 1566813 | |
| 93695 | 93695 | 231 | 0 | 0 | 783992 | -2.1429 | 0 | 0 | 1566828 | |
| 93696 | 93696 | 231 | 0 | 0 | 784001 | -2.1429 | -1 | 0 | 1566844 | |
| 93697 | 93697 | 230 | 0 | 0 | 784009 | -2.8571 | 0 | -1 | 1566860 | |
| 93698 | 93698 | 230 | 0 | 0 | 784017 | -2.8571 | 0 | 0 | 1566875 | |

- Fix the 2nd column (recover the missing events) in `A3` (the position of the car and behavioral events) in the designated row. Use the trajectory plot for help.



- Begin pre-processing

III. Pre-processing




- Reduce the amount of data
 - ◆ Take notes for preprocessing (Preprocessing log) till Part IV ends
 - Basic contents:
 - Original latency of datasets
 - # of trials before and after artifacts removal

- Components related to brain processes (with cluster ID)
- Use EEGLAB toolbox (till Part IV ends)
 - Set the path in MATLAB:



```
>> addpath(genpath('/home/eeglab5.03/'));
```
 - ```
>> eeglab
```

### III.1 Loading raw EEG data


- MATLAB function: `pop_loadcnt()` / `DR_pop_loadcnt()`
  - `DR_pop_loadcnt()`: with the same function, but doesn't read event file (actually ev2 file) when loading EEG data; use with `DR_loadcnt()`
  - EEGLAB GUI interface:
    - ◆ File → Import data → From Neuroscan .CNT file




- ◆ Switch to the desired directory and choose the CNT file (raw EEG data)




- ◆ Choose the proper resolution of raw data (currently 16-bit)



- ◆ Wait for finishing loading data...
- ◆ It takes to much time in reading the original events in raw EEG!!




- Edit name to identify datasets (there may be more than one dataset in an EEGLAB window)




- Finally...






- ◆ In command line: show the index of dataset show
- ◆ In EEGLAB GUI interface: show the index, name, and other information of the raw data
- If using EEGLAB GUI interface or `loadcnt()` → `>> EEG.event = [];`
- Remove unnecessary channels in EEG data
  - ◆ Currently (30-channel data) → remove channels with labels EKG1, EKG2, VEOU, VEOL, A1, and A2(channel indices: 1, 2, 5, 6, 27, and 33)
  - ◆ MATLAB function: `pop_select()`
  - ◆ EEGLAB GUI interface: Edit → Select data




- Type indices of channels to be REJECTED in the “Channel range” editbox, and CHECK the checkbox; otherwise, type indices of channels to be KEPT, and UNCHECK the checkbox





### III.2 Filtering and down-sampling

- Filtering: 0.5-Hz high-pass filter (HPF) and 50-Hz low-pass filter (LPF)
  - MATLAB function: `pop_eegfilt()`
  - EEGLAB GUI interface: Tools → Filter the data → Basic FIR filter




- Apply HPF and LPF SEPERATELY
  - ◆ Input the pass band (0.5) of the HPF
  - ◆ After applying the HPF, input the pass band (50) of the LPF






- ◆ Overwrite the dataset after applying each filter





If using the command-line mode, use `pop_newset()`

- Down-sampling: down-sample the sampling rate to 250 Hz
  - MATLAB function: `pop_resample()`
  - EEGLAB GUI interface:
    - ◆ Tools → Change sampling rate



- Overwrite the dataset in memory
- After this step, save the data as `sxx_yymmdd_pre.set`

- MATLAB function: `pop_saveset()`
- EEGLAB GUI interface:
  - ◆ File → Save current dataset as




### III.3 Integrating behavioral information into EEG data



- Currently, no single MATLAB code (doesn't matter)
- Procedures:
  - Load event file: `sxx_yymmdd_event.txt`
  - Extract indices of event marks (1<sup>st</sup> column in event file), event marks and trajectory (2<sup>nd</sup> column), and the latency (5<sup>th</sup> column) from the imported event file into the array `evt`
    - ◆ I. e., columns of `evt`: 1<sup>st</sup> → indices, 2<sup>nd</sup> → event marks and trajectory, 3<sup>rd</sup> → latency
  - Extract event marks (deviation onset, response onset, and response offset)

and the corresponding latencies from `evt` into the array `evt_mark`

- Save the latencies of events into `epoch_inf`. Columns:
  - 1: index of trial
  - 2~4: latencies of `dev_on`, `act_on`, and `act_off` (the PREVIOUS epoch)
  - 5~7: latencies of `dev_on`, `act_on`, and `act_off` (THIS epoch)
  - 8~10: latencies of `dev_on`, `act_on`, and `act_off` (the NEXT epoch)
  - 11: kept (1) or rejected (0); set this column to 1
 Note: latencies: the values in `evt`
- Save the RTs into `RT_original`. Columns:
  - 1: RT of this trial (sec)
  - 2: latency of the `dev_on` in this trial
- Modify the values in the 2<sup>nd</sup> column in the array `evt` (change the event marks into trajectories)
- Plot the trajectory of the car (with the event marks, indices of trials, and reaction time of each trial. Save the file into `sxx_yymmdd_raw_traj.fig`



- Integrate the event marks into EEG data
  - ◆ MATLAB function: `pop_importevent()`
  - ◆ EEGLAB GUI interface:
    - File → Import event info → From MATLAB array or ASCII file




- Append events: (checked)
- Event file or array: `[evt_mark(:, [1 2]), round(evt_mark(:, 3) / 2)]`
- Input field (column) names: `index type latency`
- Number of file header lines: 0
- Time unit (sec): 0.004
- Align event latencies to data events: NaN
- Auto adjust new events sampling rate: (checked)
- Note that the `latency` column (`round(evt_mark(:, 3) / 2)`) means 0.004 sec (sampling rate: 250 Hz); if not, time unit cannot be set into 0.004


### III.4 Importing channel locations

- This step can be done between removing unnecessary channels and extract epoch
- Procedures:
  - MATLAB function: `pop_chanedit()`
  - EEGLAB GUI interface:


◆ Edit → Channel locations



◆ Look up channel locations: use default settings (press “Ok”)




◆ Import channel locations (use “Read locations”)




◆ Choose the proper channel location file

- Experiment done before 2008 or without digitized locations of electrodes: use `30ch.xyz` (standard channel locations)
- Else, load `sxx_yymmdd.xyz` (digitized channel locations)



- ◆ Use “Opt. head center” to optimize head center (use default settings)




- ◆ Before pressing Ok in pop\_chanedit() window...
  - Use “Plot 2-D” to check the loaded channel locations
  - Check if the number of channels in EEG data is the same as that in the imported channel locations



### III.5 Extracting epochs

- Procedures:
  - Extract continuous EEG data into designated epochs (current: dev\_on, act\_on, and acf\_off)
  - MATLAB function: `pop_epoch()`
  - EEGLAB GUI interface:
    - ◆ Tools → Extract epochs






- ◆ Input event types into the text box directly or press the button in the right to select event types (dev\_on: x251, x252; act\_on: x253; act\_off: x254)
- ◆ Input the limits in epochs (dev\_on: [-1 7]; act\_on: [-2 4]; act\_off: [-2 4]) in seconds
- ◆ Rename the output dataset into `sxx_yymmdd_(epoch_type)`. Can be done when the “pop\_newset” window shows





- ◆ When the “pop\_newset” window shows, UNCHECK the “overwrite” option. Do NOT overwrite the original datasets (continuous EEG data) unless finishing extracting all kinds of epochs. (figure omitted)
- ◆ Remove baseline
  - MATLAB function: `pop_rmbase()`
  - EEGLAB GUI interface:
    - Pop-up automatically after extracting epoch, or
    - Use Tools → Remove baseline




- Range: from the beginning of each epoch (-1000 in this study) to the latency of designated events (0 in this study) in milliseconds
  - MANULLY save the dataset after finishing this step. Use the same name as the extracted dataset
- ◆ Save this dataset into “`sxx_yymmdd_(epoch_type).set`”
- Repeat the above steps until all types of epochs are extracted

### III.6 Artifact removal

- Label behavioral artifacts
  - Browse through the plotted trajectory (figure shown above)
  - Principles:
    - ◆ Reject trials with extremely short RT ( $RT < 0.3$  sec)
    - ◆ Reject trials with “overshoot” or zigzag patterns (trajectories should be flat before `dev_on` and after `act_off`)
  - Save the indices of rejected trials to `rj_behave`
- Label EEG artifacts
  - Currently, only reject epochs time-lock to `dev_on`
  - Other types of epochs: reject the same epochs as those in `dev_on`
  - Browse through the signals of EEG data
    - ◆ MATLAB function: `eegplot()`
    - ◆ EEGLAB GUI interface: Plot → Channel data (scroll)



- Set the scale to 100
- Set the time range to 5 epochs in a window




- Principles:
  - ◆ Reject channels with extreme values throughout the entire or most parts of the session (the same procedure as that in “select channels”)
  - ◆ Reject epochs with severe fluctuations across most EEG channels (browsing through the data)
- Click the epoch to be rejected
- Press “REJECT” after finishing browsing and marking the epochs to be rejected

- Press NO when confirmation window pops up (do rejection manually in this study. If pressing “Yes,” EEGLAB will reject the labeled epochs)
- `>> rj_dev_on = find(EEG.reject.rejmanul);`
- Reject artifacts
  - Merge the above indices:


```
>> rj_latency2 = union(rj_behave, rj_dev_on);
```

- Modify the 11<sup>th</sup> column in `epoch_inf` (1: kept, 0: rejected, i.e. the trials/epochs labeled above)
- Save the RT (in seconds) of the remained trial into the variable `RT` (format: 1-by-`n_of_trials` array)
- Remove epochs
  - ◆ MATLAB function: `pop_select()`
  - ◆ EEGLAB GUI interface: Edit → Select data
    - Copy the contents of `rj_latency2` (using the array editor)
    - Paste the copied data into the “Epoch range” editbox and check the checkbox



- Name the new dataset “`sxx_ymmdd_(epoch_type)_rj,`” and use the same name to save this dataset


## IV. ICA decomposition and dipole fitting



### IV.1 ICA decomposition

- Decompose EEG signals into independent signals
- Currently, only perform on epochs time-locking to `dev_on`

- MATLAB function: `pop_runica()`
- EEGLAB GUI interface: Tools → Run ICA




- Settings:
  - ◆ Running extended ICA
  - ◆ Stop criteria: maximum steps = 1024 and weight change < 1E-7




- ◆ The contents in “Command line options” editbox: (sequence doesn't matter): 'extended', 1, 'stop', 1024, 'maxsteps', 1E-7

- Run ICA
  - ◆ Check:




- ◆ Frames = (# of epochs) × (length of epochs) × (sampling rate)
- ◆ Using extended ICA
- ◆ Stop criteria: weight change < 1e-7 or after 1024 steps
- ◆ It is fine if ICA training start again
- ◆ Wait for running ICA...

- After ICA training is finished...
  - ◆ Save the dataset into “sxx\_yymmdd\_dev\_on\_rj” (the same name)
  - ◆ Copy the whole process in ICA training and save it into `sxx_yymmdd_ICA_log.txt`
  - ◆ An alternative step:
    - Before ICA training, `>> diary('sxx_yymmdd_ICA_log.txt');`
    - After ICA training, `>> diary off`
  - ◆ Plot the step vs. weight change plot (no code, using log-scale in y axis) and save it into `sxx_yymmdd_ICA_log.png`. If the trace zigzags greatly, redo noise removal on the EEG data (the convergence of training possibly by chance)




(bad ICA training)




(good ICA training)

## IV.2 ICA scalp maps


- Show the activations on different cortical regions
- Plot the ICA scalp maps
  - MATLAB function: `pop_topoplot()`
  - EEGLAB GUI interface: Plot → Component maps → In 2-D








- ◆ Set plotting electrodes to off (as above) if using standard channel locations, or set to on if using digitized channel locations




- ◆ Save the plot into “sxx\_yymmdd\_component\_map\_intrial.png”

### IV.3 Dipole fitting


- Approximate the source of EEG signals in the brain
- Currently, only perform on epochs time-locking to dev\_on
- Register channels
  - MATLAB function: `pop_dipfit_settings()`
  - EEGLAB GUI interface: Tools → Locate dipoles using DIPFIT 2.x → Head model and settings



◆ Check “No coreg” (no coregistration)




- Coarse fit: separates the brain into grids and locates dipole(s) onto appropriate grids
  - MATLAB function: `pop_dipfit_gridsearch()`
  - EEGLAB GUI interface: Tools → Locate dipoles using DIPFIT 2.x → Coarse fit (grid scan)




◆ Settings:

- Component(s): all (the speed isn't faster if just fitting some components)
- Grid in X-/Y- direction: `linspace(-85, 85, 35)` or `-85 : 5 : 85`
- Grid in Z-direction: `linspace(0, 85, 18)` or `0 : 5 : 85`  
(That is, length of each side of a grid is 5 mm)
- Rejection threshold RV (%): 40 (default)




- Fine fit: locates dipoles into appropriate grids
  - Fit only the components related to brain processes (e.g., in the above component map, only fit components 2, 3, 4, 6, 7, 8, 9, 10, 11, and 21)
  - MATLAB function: `pop_dipfit_nonlinear()`
  - EEGLAB GUI interface: Tools → Locate dipoles using DIPFIT 2.x → Coarse fit (grid scan)




◆ Basic settings:

- Component to fit: the desired component (one component/once)
- Dipole (#1 or #2): the dipole to fit. Currently, only bilateral and tangential occipital components need to fit two dipoles; otherwise, fit only one dipole
- Symmetry constrain for dipole...: check if fit two dipoles (if fitting only one dipole, it doesn't matter if this checkbox is checked)




- After fine fitting, notice the value in the “Residual variance” field should be less than the original and the value in the “position” field is usually fractions



- ◆ Use “Plot dipole(s)” to see if the dipole(s) is/are located reasonably. If

not, redo fine-fitting on this component



- Repeat the above steps unless finishing fine-fitting all the components related to brain processes
- Save the dataset into “sxx\_yymmdd\_dev\_on\_rj”

#### IV.4 Apply ICA weight matrices and dipole locations on other types of epochs

- Currently, only the dev\_on epoch performs ICA decomposition and dipole fitting
- Other types of epochs: directly apply the results of those in the dev\_on epoch
- ICA weight matrices:
  - EEG.icaact
  - EEG.icawinv
  - EEG.icasphere
  - EEG.icaweights
  - EEG.icachansind
- The dipole location matrix
  - EEG.dipfit
- Procedures
  - Load sxx\_yymmdd\_dev\_on\_rj.set. This dataset should have ICA weight matrices and dipole locations (if not, perform ICA decomposition and dipole fitting)
  - Use other variables to copy ICA weight matrices and dipole locations:

```
>> % icaact = EEG.icaact; %no need to copy this
>> icawinv = EEG.icawinv;
>> icasphere = EEG.icasphere;
```

- ```
>> icaweights = EEG.icaweights;
>> icachansind = EEG.icachansind;
>> dipfit = EEG.dipfit;
```
- **Load the dataset with other epoch:** `sxx_yymmdd_(epoch_type)_rj.set`
 - **Copy the values in the above variables**

```
>> % EEG.icaact = icaact;   %%%WRONG!!
>> EEG.icaact = icaweights * icasphere * ...
    reshape(EEG.data, size(EEG.data, 1), ...
    size(EEG.data, 2) * size(EEG.data, 3));
>> EEG.icaact = reshape(EEG.icaact, ...
    size(EEG.data, 1), size(EEG.data, 2), ...
    size(EEG.data, 3));
>> EEG.icawinv = icawinv;
>> EEG.icasphere = icasphere;
>> EEG.icaweights = icaweights;
>> EEG.icachansind = icachansind;
>> EEG.dipfit = dipfit;
```
 - `>> eeglab redraw;` %upgrade the information in the dataset
 - **Save the dataset into** `sxx_yymmdd_(epoch_type)_rj.set`
 - Repeat the above steps until all types of epochs have ICA weight matrices and dipole locations

V. Automatically preprocessing

V.1 Why using MATLAB scripts for pre-processing

- Drawbacks of using EEGLAB menu to pre-process EEG data
 - Time-consuming
 - ◆ Wait for finishing one step
 - ◆ Wait for user's response
 - Make mistakes easily
 - ◆ Type wrong numbers, press wrong key, ...
 - Arduous and inefficient
- Since you won't pre-process/analysis only one EEG dataset, why not using/writing MATLAB scripts?
- Currently, MATLAB script for pre-process is available


V.2 MATLAB scripts for pre-processing

- Script: `DR_preprocessing_compact.m` (ver 200809)
 - Pre-process EEG data semi-automatically

- ◆ With several breakpoints to prevent from overwriting data
- ◆ Some steps still need manual operations
- Can process multiple datasets
- Integrate pre-process of behavioral and EEG data in a whole (the parts III and IV in this document)
- Required toolbox/scripts:
 - ◆ EEGLAB ver. 5.03
 - ◆ `mypath.m` (set the path in MATLAB)
 - ◆ `DR_pop_loadcnt.m`
 - ◆ `DR_loadcnt().m`
 - ◆ Still have to make preprocessing log (contents described before)
- Check before you run:
 - ◆ The path is correct or not
 - ◆ Select or add dataset in the variable `SET`
 - ◆ Check if the output directory has other files. Files with the same names as the output files will be **OVERWRITTEN!**
- The following steps need manual operation: (so check the progress occasionally)
 - ◆ Artifacts removal (both behavioral and EEG data)
 - ◆ Plot ICA weight change
 - ◆ Fine-fitting dipoles
- Outputs
 - ◆ Datasets
 - `sxx_yymmdd_pre.set`: filtered and down-sampled EEG data
 - `sxx_yymmdd_(epoch_type).set`: epoched EEG data
 - `sxx_yymmdd_(epoch_type)_rj.set`: epoched and artifact-removed EEG data with ICA weights and dipole locations
 - ◆ MAT files
 - `sxx_yymmdd_epoch_inf.mat`
`epoch_inf`: time frames of the events in a trial, the previous trial, and the next trial, and a flag to indicate if the trial is kept or not
 - `sxx_yymmdd_RT_of_trials.mat`
`RT`: reaction time (sec) of the remained trials after artifact removal
`RT_original`: reaction time (sec) of the original trials
 - `sxx_yymmdd_rejected_trials.mat`
`rj_behave`: trials behavioral artifacts
`rj_dev_on`: epochs with artifacts in EEG data
`rj_latency2`: union of the above two variables
`rj_latency1`: (obsolete, kept for future use)
 - ◆ Figures (shown in part III and IV)


- `sxx_yymmdd_raw_traj.fig`: raw trajectory in a session
- `sxx_yymmdd_ICA_log.png`: log file for ICA training [plot manually]
- `sxx_yymmdd_component_map_intrial.png`: component map after ICA decomposition
- ◆ Txt file
 - `sxx_yymmdd_preprocess_log.txt`: preprocessing log file [create and write the contents manually]
 - `sxx_yymmdd_ICA_log.txt`: log file for ICA training
- Knowing bugs
 - ◆ Bugs in applying ICA weight matrices and dipole locations on other types of epochs. The epochs may not have the desired ICA weight matrices and dipole location. Check after pre-process ends

VI. Computation of tonic power spectra



- Script: `DR_PBase_setup.m`
 - Save the information about components from all datasets and settings for plotting figures
 - ◆ The “database” of subjects and components (for checking)
 - Add or remove datasets/components if needed
 - ◆ Be cautious!
 - ◆ Components and clusters #:
 - Clusters #: use the settings in `DR_PBase_setup.m`
 - Pre-classify when pre-process ends
 - Modify after clustering (described in Section VII.2): make the cluster # negative (e.g.: 5 → -5)

- All processes below need this script
- Script: `RS_powerbase.m` (ver 200901)
 - Compute and plot power spectra for SINGLE-SUBJECT data
 - Two parts in this script:
 - ◆ Compute power spectra
 - Flowchart:



- ◆ Plot power spectra
- Required scripts
 - ◆ `mypath.m`
 - ◆ `DR_PBase_setup.m`
 - ◆ `timefreq2004.m`: performs time-frequency transform. The same as `timefreq.m` in UCSD, rename to prevent confusing in MATLAB
- Required inputs
 - ◆ Components and identified clusters (in `DR_PBase_setup.m`)
 - ◆ ICA activation of the above components (in EEG dataset)
- Check before you run:
 - ◆ Variables
 - `PLOT_SET`: the plotted datasets
 - `MN: conditions.` ('motionless' or 'motion')
 - `p_val`: significance level (currently $1E-4$)
 - `epoch_type`: types of power spectra
 - ' ' → tonic spectra
 - 'dev_on' → phasic spectra
 - 'all' → mixed spectra
 - `rj: '_rj'` → artifacts removed dataset
 - ◆ Directories
 - Tonic power spectra: `FilePath/MN/IC00~IC10`
 - Other power spectra: `FilePath/epoch_type/MN/IC00~IC10`
- Detailed process of computing power spectra
 - ◆ Perform time-frequency transform using `timefreq2004.m`
 - `[PB, freqs, times] = timefreq2004(baseline, SR, 'wavelet', 0, 'freqs', [min_freq max_freq], 'winsize', 128, 'pdratio', 2, 'ntimesout', 100);`

baseline: ICA activations in (part of) an epoch; not necessarily the activations before deviation

- ◆ Compute magnitude and convert into dB
 - `PB = 10 * log10(PB .* conj(PB));`
- ◆ Average the resulting spectra in different time windows and convert into a 2-D array
 - `PB_mean = trimmean(PB, 10, 2);`
 - `PB_mean = reshape(PB_mean, size(freqs, 2), n_of_trials);`
- ◆ Choose desired frequency range
 - `PB_mean = PB_mean(find(freqs >= 3 & freqs <= 45), valid_trial) %valid_trial: trials with desired RTs`
 - `freqs = freqs(find(freqs >= 3 & freqs <= 45));`
- ◆ Sort by RT
 - `PB_mean = PB_mean(:, ur_idx);`
- ◆ Compute alert baseline spectra (spectra from trials with shortest 10% RT)
 - `PB_alert = PB_mean(:, 1 : alert.trials);`
 - `PB_alert_mean = trimmean(PB_alert, 10, 2) * ones(1, size(valid_trial, 2));`
- ◆ Normalize
 - `PB_n = PB_mean - PB_alert_mean;`
- ◆ Moving average (RT and power spectra, using for loop)
 - Power spectra: `PB_mov = [PB_mov trimmean(PB_n(:, k : k + alert.trials - 1), 10, 2)];`
 - RT: `RT_s_mov = [RT_s_mov trimmean(RT_s(:, k : k + alert.trials - 1), 10, 2)];`
- Detailed process of applying two-sampled *t*-test (using for loop)
 - ◆ A moving window
 - `PB_tmp = PB_mean(:, sig_idx(k) : sig_idx(k) + alert.trials - 1);`
 - ◆ Estimate significance level
 - `[H, P] = ttest2(PB_tmp', PB_alert_mean(:, 1 : alert.trials)', p_val / 2 / size(freqs, 2), 'both');`
 - `PB_alert_mean(:, 1 : alert.trials):` alert baseline spectra
 - `p_val / 2 / size(freqs, 2):` *p* value corrected by degree of freedom
 - H: show if the power on this frequency bin is significant
 - P: *p* value on this frequency bin
 - ◆ Merge into the matrix (for adding contour on power image)


- `H_mask(:, k) = H';`
- `P_mask(:, k) = P';`

■ Outputs


◆ MAT file

- `sxx_yymmdd_PBaseRT_(component #)(rj).mat`
`icawinv`: ICA inverse weight matrix of this component
`chanlocs`: channel locations of this component
`PB`: original power spectra computed by `timefreq2004.m`
`PB_mean`: power spectra (averaged across time windows)
`PB_alert`: power spectra of alert trials
`PB_alert_mean`: averaged power spectra of alert trials (alert baseline spectra)
`PB_n`: normalized power spectra
`PB_mov`: moving averaged power spectra (from `PB_n`)
`freqs`: output frequency bins
`RT`: reaction time of all trials
`RT_s`: sorted reaction time
`RT_s_mov`: moving averaged reaction time
`FREQ_INC`: power spectra of the extracted frequency bands
`H_mask`: significant region of power spectra
`P_mask`: *p* value of individual frequency/trial
`p_val`: the set *p* value (uncorrected)
`URPB`: (sparse moving average)
`URPB_mov`: (sparse moving average)


- ◆ Figures: [if plotting tonic power spectra, this script plots in `fig` and `png` files; else, only plots in `png` format]
 - Power spectra after each step: `sxx_yymmdd_PBaseRT_(component #)(rj)_method.png/fig`



- Power image and trends of power increase on the frequency bands defined in `DR_PBase_setup.m: sxx_yymmdd_PBaseRT_(component #) (rj)_result.png/fig`




- The same plots in the above point and sparse moving average on the previous bands: `sxx_yymmdd_PBaseRT_(component #) (rj)_sparse.png/fig`



- Knowing bugs: (not reported)


VII. Clustering and component re-selection



- Clustering: characterize the common pattern of spectral activities of similar ICs from different subjects and sessions
- Component re-selection: remove abnormal spectra/component scalp maps from IC clusters to refine the trend
- Re-compute cluster spectra after re-selecting components

VII.1 Clustered RTs

- Script: DR_PlotClustRT.m (ver 200903)
 - Plot CLUSTERED RTs under motionless or motion condition
 - Required scripts
 - ◆ mypath.m
 - ◆ DR_PBase_setup.m
 - Required inputs
 - ◆ RTs of trials after artifact removal
 - Check before you run:
 - ◆ Variables
 - PLOT_SET: the plotted datasets
 - MN: conditions. ('motionless' or 'motion')
 - MN_COLOR: color code used to indicate motionless or motion conditions
 - rj: '_rj' → artifacts removed dataset
 - ◆ Directories: FilePath/MN/
 - Outputs
 - ◆ MAT file: RT_Sorted_Trial_(MN).mat
 - urRT_s: sorted clustered RTs (before artifacts removal)
 - RT_s: sorted clustered RTs (after artifacts removal)
 - MN: motionless or motion condition
 - session_count: # of sessions in this cluster
 - ◆ Figures: RT_Sorted_Trial_(MN).fig. Print the png file in PC



VII.2 Clustering ICA scalp maps, alert baseline power, and dipole lo-

cations

- **Script:** DR_PlotClustCompPowerDipole.m (ver 200903)
 - **Compute and plot averaged scalp map, alert baseline power, and dipole locations in a cluster**
 - **Required scripts**
 - ◆ mypath.m
 - ◆ DR_PBase_setup.m
 - **Required inputs**
 - ◆ ICA activation of the components in EEG datasets
 - ◆ Identified clusters with components from EEG datasets
 - ◆ Standard channel locations (for finding missing channels in datasets)
 - ◆ Power spectra of the components in the plotted cluster
 - ◆ Dipole locations of the components in the plotted cluster (optional)
 - **Check before you run:**
 - ◆ **Variables**
 - PLOT_SET: the plotted datasets
 - MN: conditions. ('motionless' or 'motion')
 - MN_COLOR: color code used to indicate motionless or motion conditions
 - n_of_cls: # of cluster; undef_cls: undefined cluster
 - rj: '_rj' → artifacts removed dataset
 - std_chanlocs: file of standard channel locations (for checking missing channels)
 - avg: how to compute average. '_trim': using 10% trimmed mean
 - ◆ **Directories:** FilePath/MN/
 - **Detailed process of computing average scalp map in a cluster**
 - ◆ Automatically negate the weights if needed to make all the clusters in this cluster with similar distribution in color: normalize each weight column vector according to it maximum absolute value
 - Determine if each scalp map in this cluster has missing channels

```
[tmp_chanlabel chanlabel_diff] =  
setdiff(upper(std_chanlabel), upper(chanlabel));
```
 - If any dataset in this cluster has missing channels, fill the corresponding positions in ICA inverse weight vector with NaN, and then negate the weight if needed

```
icawinv_tmp(find(~isnan(icawinv_tmp))) =  
EEG.icawinv(:, plot_comp);  
[abs_max max_idx] = max(abs(icawinv_tmp));  
icawinv      =      [icawinv      icawinv_tmp      /
```

```
icawinv_tmp(max_idx)];
```

- If with no missing channels, negate the weights if needed

```
[abs_max max_idx] =
```


```
max(abs(EEG.icawinv(:, plot_comp)));
```

```
icawinv = [icawinv (EEG.icawinv(:, plot_comp)) /  
(EEG.icawinv(max_idx, plot_comp))];
```

- ◆ Collect inverse weight matrix in a cluster

- `chan all = [chan all chanlocs];`

- ◆ Plot each scalp map in this cluster for manually negating the weights (enter the indices in the MATLAB command line in the matrix form)



- ◆ Compute average inverse weight vector in this cluster

- If any dataset in this cluster has missing channels, compute the average vector row by row

```
icawinv tmp = find(~isnan(icawinv(k, :)));
```

```
icawinv tmp = mean(icawinv(k, icawinv tmp));
```

```
avg icawinv = [avg icawinv; icawinv tmp];
```

- If with no missing channels, compute the average vector directly

```
avg icawinv = mean(icawinv, 2);
```

- ◆ Compute average paired correlation in this cluster

- Compute paired correlation and take upper triangle matrix of the resulting matrix


```
comp_corr = triu(corrcoef(icawinv));
```

- **Set the diagonal value to 0 (original 1)**

```
diag_matrix = zeros(size(comp_corr));
for k = 1 : size(comp_corr, 1)
    diag_matrix(k, k) = comp_corr(k, k);
end
```

- **Compute average and SD of correlation value**

```
tmp = find(comp_corr > 0);
avg_corr = mean(comp_corr(tmp));
sd_corr = std(comp_corr(tmp));
```

■ Detailed process of computing average dipole locations

- ◆ **Fit dipoles location if the input dataset does not have dipole locations (doing when loading datasets; the same procedures as those in the pre-processing part)**

```
● if ~isfield(EEG, 'dipfit') || isempty(EEG.dipfit)
    .....
end
```

- ◆ **Using dipplot() to compute the Talairach coordinates of each dipole; save in dip_tmp for further use**

```
● eval(...
    ['dip_tmp = dipplot(EEG.dipfit.model(plot_comp),
    'mri', mrifile, 'view', view_angle, 'color',
    {[1 1 0]}], ' ...
    ''dipolesize'', 20, ''dipolelength'', 0, '
    [dip_option ...
    ', 'projimg'', 'off'', 'projlines'', 'on'',
    'coordformat'', 'spherical'''] ');]);
```

- ◆ **Update and save the original EEG dataset**

```
● EEG.dipfit.model(plot_comp).eleccoord
    = dip_tmp.eleccoord;
● EEG.dipfit.model(plot_comp).mnicoord
    = dip_tmp.mnicoord;
● EEG.dipfit.model(plot_comp).talcoord
    = dip_tmp.talcoord;
● EEG = pop_saveset(EEG, [FilePath plot_set SL
    plot_set EpochType{1} rj '.set']);
```

- ◆ **Collect dipole locations and compute total residual variance in a cluster**

```
● dipole_all = [dipole_all dip_tmp];
● rv_all = rv_all + dip_tmp.rv;
```

- ◆ **Collect Talairach coordinates and other information in a cluster [use**

`talcoord_all` (Talairach coordinates of dipoles) for example]

- If the cluster is the occipital or tangential occipital cluster (possibly with two dipoles)

```
if dip_tmp.talcoord(k, 1) < 0 %left
    talcoord_all(j, 1:3) = dip_tmp.talcoord(k, :);
    .....
else %right
    talcoord_all(j, 4:6) = dip_tmp.talcoord(k, :);
    .....
end
```

- If the cluster is one of the other clusters (only one dipole)

```
talcoord_all(j, :) = dip_tmp.talcoord(1, :);
```

- ◆ Compute mean and SD of residual variance and Talairach coordinates in a cluster [use `talcoord_all` for example]

- `rv_all = rv_all / cls_length; %rv`
- `talcoord_all(j + 1, k) = mean(talcoord_all(find(~isnan(talcoord_all(1 : j, k))), k)); %average`
- `talcoord_all(j + 2, k) = std(talcoord_all(find(~isnan(talcoord_all(1 : j, k))), k)); %SD`

- `diffmap`: unknown function, keep in the script

- ◆ Knack for plotting dipole locations in the same axis

- The first dipole (or pair of dipoles)
Plot dipoles on the designated axis directly
Then, set this axis as the target (for pasting) (`target = gca;`)
- The other (pairs of) dipoles:
Plot dipoles on ANOTHER axis (dummy figure)
Then, find and copy the graphics objects with “line” properties in this axis/figure (`h2 = findobj(gcf, 'Type', 'line');`)
Next, paste these objects on the designated axis (`copyobj(h2, target);`)
Finally, close the dummy figure

- Detailed process of computing alert baseline spectra

- ◆ Take transposition of alert baseline spectra of each input dataset in this cluster to facilitate plotting spectra traces


- `PB_alert_mean = PB_alert_mean(:, 1)';`

- ◆ Collect the spectra

- `spec_all = [spec_all; PB_alert_mean];`


- Outputs

- ◆ **MAT file:** `ICnn_component_(cls_ver)_(MN).mat` (nn: code of this cluster)
 - `std_chanlocs`: standard channel locations
 - `chan_all`: channel locations of all datasets in this cluster
 - `icawinv`: ICA inverse weight matrix of all datasets in this cluster
 - `avg_icawinv`: average ICA inverse weight matrix of all datasets in this cluster
 - `comp_corr`: paired correlation of the components in this cluster
 - `avg_corr`: average correlation of the components in this cluster
 - `sd_corr`: SD of correlation of the components in this cluster
 - `spec_all`: alert baseline spectra from the datasets in this cluster
 - `freqs`: frequency bins
 - `dipole_all`: information of all dipoles in this cluster
 - `dipole_avg`: information of average dipoles in this cluster
- ◆ **txt file:** `ICnn_ClsDipole_(cls_ver)_(MN).txt`
 - Save the Talairach coordinates in this cluster in a printable file
- ◆ **Figures:**
 - Average scalp maps in this cluster:




- Alert baseline spectra in this cluster:

`ICnn_Cluster_Spectra_(cls_ver)_(MN).fig/png`




- Dipole locations in this cluster: (in array of plots)

ICnn_DipoleMap_(cls_ver)_(MN).fig/png



- Dipole locations in this cluster: (in a figure)

ICnn_DipoleMap_(cls_ver)_single_MN.fig/png



- Knowing bugs:
 - ◆ Having errors when running this script in other accounts. Checking
- Script: DR_PlotClustCompmap.m (ver 200904)
 - Plot average component map in clusters
 - Read ICA inverse weight vectors and plot figures
 - Required scripts
 - ◆ mypath.m
 - ◆ DR_PBase_setup.m
 - Required inputs
 - ◆ ICA activation of the above components (in EEG dataset)
 - ◆ Components and identified clusters
 - Check before you run:
 - ◆ Variables
 - MN: conditions
 - MN_COLOR: color code used to indicate motionless or motion conditions
 - p_val: significance level (currently 1E-10)
 - rj: '_rj' → artifacts removed dataset
 - ◆ Directories
 - Tonic power spectra: FilePath/MN/
 - Other power spectra: FilePath/epoch_type/MN/
 - Outputs
 - ◆ Figure: Avg_Scalp_Map_(cls_ver)_(MN).fig. Print the png file in PC



- Knowing bugs: (not reported)

VII.3 Clustering power spectra

- Script: `DR_PlotClustPbase.m` (ver 200903)
 - Compute and plot power spectra for CLUSTERED data
 - The process and structure are similar to those in `RS_powerbase.m` except estimating the significance level
 - Required scripts
 - ◆ `mypath.m`
 - ◆ `DR_PBase_setup.m`
 - Required inputs
 - ◆ Components and identified clusters
 - ◆ Average ICA inverse weight vector in this cluster
 - ◆ Power spectra of the components in the plotted cluster
 - Check before you run:
 - ◆ Variables
 - `PLOT_SET`: the plotted datasets
 - `MN: conditions.` ('motionless' or 'motion')
 - `MN_COLOR`: color code used to indicate motionless or motion conditions
 - `p_val`: significance level (currently $1E-10$)
 - `epoch_type`
 - ' ' → tonic spectra
 - 'dev_on' → phasic spectra

- 'all' → mixed spectra
 - rj: '_rj' → artifacts removed dataset
 - SD:
 - '' → plot avg. of spectra
 - '_SD' → plot (avg. + SD)
 - no_band: special purpose for thesis
 - '' → plot all the bands defined in DR_PBase_setup.m
 - '1' → plot all except the 3rd band defined in DR_PBase_setup.m
 - nor: how to normalize spectra. Currently 1
- ◆ Directories
 - Tonic power spectra: FilePath/MN/
 - Other power spectra: FilePath/epoch_type/MN/
- Detailed process of normalizing power spectra
 - ◆ Collect alert baseline power and normalized power spectra from the datasets in this cluster
 - PB_alert: padded to keep the same dimension as that of PB_n


```
PB_alert = [PB_alert tmp_set.PB_alert_mean];
PB_n = [PB_n tmp_set.PB_n];
```
 - PB_alert_subj: one dataset, one spectrum


```
if ~strcmp(old_set, plot_set)
    PB_alert_subj = ...
    [PB_alert_subj tmp_set.PB_alert_mean(:, 1)];
end
```
 - ◆ After finishing loading datasets, sort baseline power by RT
 - ```
[RT_s ur_idx] = sort(RT_s, 2);
PB_mean = PB_mean(:, ur_idx);
PB_n = PB_n(:, ur_idx);
PB_alert = PB_alert(:, ur_idx);
subj_seq = subj_seq(:, ur_idx);
```
  - ◆ Moving average (power spectra and RT)
    - ```
PB_mov = [PB_mov trimmean(PB_n(:, j : j + alert.trials - 1), 10, 2)];
```
 - ```
RT_s_mov = [RT_s_mov trimmean(RT_s(:, j : j + alert.trials - 1), 10, 2)];
```
- Detailed process of applying two-sampled *t*-test (using for loop)
  - ◆ A moving window
    - ```
PB_tmp = PB_mean(:, sig_idx(k) : sig_idx(k) + alert.trials - 1);
```
 - ```
subj_seq_id_tmp = subj_seq_id(:, sig_idx(k) : sig_idx(k) + alert.trials - 1);
```

- ◆ Compute weighted mean alert baseline power in this cluster (using for loop)
  - Find # of trials from each session in this cluster  
`n_of_trial = length(find(subj_seq_id_tmp == m));`
  - If there are trials in this component (session) in this cluster, compute the weighted sum  
`PB_alert_tmp = PB_alert_tmp + PB_alert_subj(:, m) * n_of_trial;`
  - Update # of subjects (if this dataset didn't appear)  
`n_of_subj_win = n_of_subj_win + 1;`
  - Compute weighted average of alert baseline spectra  
`PB_alert_tmp = PB_alert_tmp * ones(1, alert.trials) / alert.trials;`
- ◆ Estimate significance level
  - `[H, P] = ttest2(PB_tmp', PB_alert_tmp', p_val / 2 / size(freqs, 2) / n_of_subj_win, 'both');`
  - `p_val / 2 / size(freqs, 2) / n_of_subj_win`: *p* value corrected by degree of freedom
  - H: show if the power on this frequency bin is significant
  - P: *p* value on this frequency bin
- ◆ Merge into the matrix (for adding contour on power image)
  - `H_mask(:, k) = H';`
  - `P_mask(:, k) = P';`
- Outputs
  - ◆ MAT file:
    - Tonic spectra:  
`ICnn_PBase_(nor)_(MN).mat`
    - `cls_length`: # of components in this cluster
    - `session_count`: # of sessions in this cluster
    - `icawinv`: averaged ICA inverse weight matrix of this cluster
    - `chanlocs`: averaged channel locations of this component
    - `freqs`: output frequency bins
    - `RT_s`: sorted reaction time in this cluster
    - `RT_s_mov`: moving averaged reaction time in this cluster
    - `subj_seq`: sequence of subject (corresponding to each trial)
    - `p_val`: the set *p* value (uncorrected)
    - `PB_mean`: power spectra (averaged across time windows)
    - `PB_alert`: averaged power spectra of alert trials
    - `PB_n`: normalized power spectra
    - `PB_alert_subj`: power spectra of alert trials (from the original



dataset of each subject)

PB\_mov: moving averaged power spectra (from PB\_n)

FREQ\_INC: power spectra of the extracted frequency bands

H\_mask: significant region of power spectra

P\_mask:  $p$  value of individual frequency/trial

IC (nn) : datasets and components

- Other spectra:

ICnn\_PBase\_(nor)\_(MN)\_(epoch\_type).mat

(variables are the same as those in ICnn\_PBase\_nor\_MN.mat)


- ◆ Figures: [if plotting tonic power spectra, this script plots in fig and png formats. (Note: print the png file again in PC); else, only plots in png format]

- Power spectra after each step:

ICnn\_PBase\_(nor)\_(MN)\_method.png/fig for tonic spectra

ICnn\_PBase\_(nor)\_(MN)\_(epoch\_type)\_method.png/fig

for other types of spectra



- ICnn\_PBase\_(nor)\_(MN)\_result.png/fig or

ICnn\_PBase\_(nor)\_(MN)\_(epoch\_type)\_result.png/fig:

(plotted by DR\_PlotClustPbase\_merge.m, described later)


- Knowing bugs: (not reported)

- Script: DR\_PlotClustPbase\_merge.m (ver 200903)

- Plot power spectra for CLUSTERED data with dipole locations

- Run after finishing plotting clustered dipole locations (using script DR\_PlotClustCompPowerDipole.m, described in Section VII.2)

- Required scripts
  - ◆ mypath.m
  - ◆ DR\_PBase\_setup.m
- Required inputs
  - ◆ Average ICA inverse weight vector in this cluster
  - ◆ Clustered power spectra of this cluster
  - ◆ Dipole locations and alert baseline spectra of this cluster
- Check before you run:
  - ◆ Variables
    - PLOT\_SET: the plotted datasets
    - MN: conditions. ('motionless' or 'motion')
    - MN\_COLOR: color code used to indicate motionless or motion conditions
    - p\_val: significance level (currently 1E-10)
    - epoch\_type: types of power spectra
      - ' ' → tonic spectra
      - 'dev\_on' → phasic spectra
      - 'all' → mixed spectra
    - rj: '\_rj' → artifacts removed dataset
    - SD: ' ' → plot avg. of spectra; '\_SD' → plot (avg. + SD)
    - no\_band: special purpose for thesis
      - ' ' → plot all the bands defined in DR\_PBase\_setup.m
      - '1' → plot all except the 3<sup>rd</sup> band defined in DR\_PBase\_setup.m
    - nor: how to normalize spectra. Currently 1
  - ◆ Directories
    - Tonic power spectra: FilePath/MN/
    - Other power spectra: FilePath/epoch\_type/MN/
- Outputs
  - ◆ Figure: [if plotting tonic power spectra, this script plots in fig and png files (Note: print the png file again in PC); else, only plots in png format]
    - Power spectra for the results
      - ICnn\_PBase\_(nor)\_(MN)\_method.png/fig for tonic spectra [and then print the png file in PC], or
      - ICnn\_PBase\_(nor)\_(MN)\_(epoch\_type)\_result.png/fig for other types of spectra




- Knowing bugs: (not reported)

#### VII.4 Component re-selection

- Criteria
  - ICA scalp maps in this cluster
  - Alert baseline spectra in this cluster
  - Dipole locations of independent components in this cluster (not necessary)
- Modify the cluster # negative in `DR_PBase_setup.m` (e.g.: 5 → -5)
- Re-evaluate clustered RT, power spectra, alert baseline spectra, average scalp map, and dipole locations

#### VIII. Further comparisons




- Further comparisons in EEG (compare different conditions, etc.)
- The logic of scripts in this part mostly base on that of the scripts in the previous part, or just read the `MAT` files generated in the previous section and plot figures

### VIII.1 Compare clustered RTs between motionless and motion conditions

- Script: `DR_PlotClustRT_motionless_motion.m` (ver 200903)
  - Compare CLUSTERED RTs under motionless and motion condition
  - Required scripts
    - ◆ `mypath.m`
    - ◆ `DR_PBase_setup.m`
  - Required inputs
    - ◆ RTs of trials after artifact removal
  - Check before you run:
    - ◆ Variables
      - `PLOT_SET`: the plotted datasets
      - `MN`: conditions (currently `MN = {'motionless', 'motion'};`)
      - `rj`: `'_rj'` → artifacts removed dataset
      - `p_val`: significance level (currently `1E-4`; however, this is meaningless due to the appropriate method for estimating the significance level is yet to be found)
    - ◆ Directories: `FilePath/MN/`
  - Outputs
    - ◆ `MAT` file: `RT_Sorted_Trial_motionless+motion.mat`
      - `PLOT_SET1`: list of motionless datasets
      - `PLOT_SET2`: list of motion datasets
      - `urRT_s1`: sorted clustered RTs of motionless datasets (before artifact removal)

- `urRT_s2`: sorted clustered RTs of motion datasets (before artifact removal)
- `RT_s1`: sorted clustered RTs of motionless datasets (after artifact removal)
- `RT_s2`: sorted clustered RTs of motion datasets (after artifact removal)
- `RT_s`: clustered sorted reaction time (after artifacts removal)
- `RT_s_mov`: moving averaged reaction time in this cluster
- `MN`: motionless or motion condition
- `p_val`: the set  $p$  value (currently  $1E-4$ ) [meaningless]
- `norm_distr`: check if the clustered datasets are with normal distribution (the 1<sup>st</sup> element: motionless datasets; 2<sup>nd</sup>: motion datasets)
- `stat_diff`: struct for storing actual  $p$  value and flags to show if two conditions are statistically significant [meaningless]  
 $H$ : flags to show if two conditions are statistically significant  
 $P$ : actual  $p$  value

- ◆ **Figures:** `RT_Sorted_Trial_motionless+motion.fig`. Print the png file in PC



## VIII.2 Compare power spectra across motionless and motion conditions

- Script: `DR_PlotClustPbase_motionless_motion.m` (ver 200904)
  - Plot power spectra for CLUSTERED data across motionless and motion conditions

- The logic is similar to that in `DR_PlotClustPbase.m`
- Since some of the subjects may participated only motionless or motion session, it is necessary to exclude such datasets and re-cluster the data
- Required scripts
  - ◆ `mypath.m`
  - ◆ `DR_PBase_setup.m`
- Required inputs
  - ◆ Components and identified clusters
  - ◆ Average ICA inverse weight vector in this cluster **[WRONG!! Since subjects may be different, it is necessary to re-compute average ICA inverse weight matrix]**
  - ◆ Power spectra of the components in the plotted cluster
- Check before you run:
  - ◆ Variables
    - `PLOT_SET`: the plotted datasets
    - `MN`: conditions (currently `MN = {'motionless'}, {'motion'};`)
    - `MN_COLOR`: color code used to indicate motionless or motion conditions
    - `p_val`: significance level (currently `1E-10`)
    - `epoch_type`: types of power spectra
      - ' ' → tonic spectra
      - 'dev\_on' → phasic spectra
      - 'all' → mixed spectra
    - `rj`: '\_rj' → artifacts removed dataset
    - `SD`: ' ' → plot avg. of spectra; '\_SD' → plot (avg. + SD)
    - `no_band`: special purpose for thesis
      - ' ' → plot all the bands defined in `DR_PBase_setup.m`
      - '1' → plot all except the 3<sup>rd</sup> band defined in `DR_PBase_setup.m`
    - `nor`: how to normalize spectra. Currently 1
  - ◆ Directories
    - Tonic power spectra: `FilePath/MN/`
    - Other power spectra: `FilePath/epoch_type/MN/`
- Outputs
  - ◆ MAT file:
    - Tonic spectra: (nn: code of this cluster)  
`ICnn_PBase_nor_motionless+motion.mat`  
`motionless`: struct for storing necessary variables (motionless)  
`IC(nn)`: datasets and components  
`session_count`: # of sessions in this cluster  
`cls_length`: # of components in this cluster

`icawinv`: averaged ICA inverse weight matrix of this cluster  
`chanlocs`: averaged channel locations of this cluster  
`freqs`: output frequency bins  
`RT_s`: sorted reaction time in this cluster  
`RT_s_mov`: moving averaged reaction time in this cluster  
`RT_s_overlapidx`: indices of trials with overlapped RT  
`RT_s_mov_overlapidx`: indices of trials with overlapped moving averaged RT  
`subj_seq`: sequence of subject (corresponding to each trial)  
`PB_mean`: power spectra (averaged across time windows)  
`PB_alert`: averaged power spectra of alert trials  
`PB_alert_subj`: power spectra of alert trials (from the original dataset of each subject)  
`PB_n`: normalized power spectra  
`PB_mov`: moving averaged power spectra (from `PB_n`)  
`FREQ_INC`: power spectra of the extracted frequency bands  
`H_mask`: significant region of power spectra  
`P_mask`: *p* value of individual frequency/trial  
`motion`: (the same function as that in the above one)  
`RT_s_overlap`: overlapped RT in both conditions  
`RT_s_mov_overlap`: moving averaged overlapped RT in both conditions

- Other spectra:


`ICnn_PBase_(nor)_motionless+motion_(epoch_type).mat`

(variables are the same as those in `ICnn_PBase_nor_MN.mat`)

- ◆ Figures: if plotting tonic power spectra, this script plots in `fig` and `png` files [Note: print the `png` file again in PC]; else, only plots in `png` format

- Power spectra after each step:

`ICnn_PBase_(nor)_motionless+motion.png/fig` (tonic),  
`ICnn_PBase_(nor)_motionless+motion_(epoch_type)_method.png/fig` (other)



- Knowing bugs: (not reported)

### VIII.3 Compare tonic and mixed power spectra

- Script: DR\_PlotClustPBase\_compare.m (ver 200905)
  - Compare power spectra between Two of the states: tonic, mixed, and phasic
  - Read the power spectra and plot figures
  - Required scripts
    - ◆ mypath.m
    - ◆ DR\_PBase\_setup.m
  - Required inputs
    - ◆ Components and identified clusters
    - ◆ Average ICA inverse weight vector in this cluster
    - ◆ Power spectra of the plotted cluster
  - Check before you run:
    - ◆ Variables
      - PLOT\_SET: the plotted datasets
      - MN: conditions
      - MN\_COLOR: color code used to indicate motionless or motion conditions
      - p\_val: significance level (currently 1E-10)
      - EPOCH\_TYPE: types of power spectra (2 of the followings; currently EPOCH\_TYPE = {'', 'all'};)
      - '' → tonic spectra
      - 'dev\_on' → phasic spectra



'all' → mixed spectra

- EPOCH\_COLOR: color code used to indicate different types of spectra
- rj: '\_rj' → artifacts removed dataset
- SD: '' → plot avg. of spectra; '\_SD' → plot (avg. + SD)
- no\_band: special purpose for thesis
  - '' → plot all the bands defined in DR\_PBase\_setup.m
  - '1' → plot all except the 3<sup>rd</sup> band defined in DR\_PBase\_setup.m
- nor: how to normalize spectra. Currently 1


#### ◆ Directories

- Tonic power spectra: FilePath/MN/
- Other power spectra: FilePath/epoch\_type/MN/

#### ■ Outputs

##### ◆ Figure:

ICnn\_compare\_(epoch\_type1)\_(epoch\_type2)\_(MN) (nor)  
(rj) (SD) (noband).fig. Print the png file in PC



#### ■ Knowing bugs: (not reported)

- Script: DR\_PlotClustPBase\_powinc.m (ver 200906) [only use in defense slide]
- Compare power increase between
  - ◆ Two of the following states: tonic, mixed, and phasic,
  - ◆ Motionless and Motion conditions, and
  - ◆ Different clusters
- Plot power on a frequency band in a figure

- Read the power spectra and plot figures
- Required scripts
  - ◆ mypath.m
  - ◆ DR\_PBase\_setup.m
- Required inputs
  - ◆ Components and identified clusters
  - ◆ Average ICA inverse weight vector in this cluster
  - ◆ Power spectra of the plotted cluster
- Check before you run:
  - ◆ Variables
    - PLOT\_SET: the plotted datasets
    - MN: conditions. (currently MN = {'motionless', 'motion'};)
    - p\_val: significance level (currently 1E-10)
    - EPOCH\_TYPE: types of power spectra (2 of the followings; currently EPOCH\_TYPE = {'', 'all'};)
    - '' → tonic spectra
    - 'dev\_on' → phasic spectra
    - 'all' → mixed spectra
    - rj: '\_rj' → artifacts removed dataset
    - SD: '' → plot avg. of spectra; '\_SD' → plot (avg. + SD)
    - no\_band: special purpose for thesis (not important in this script since each figure only has power in a frequency band)
    - '' → plot all the bands defined in DR\_PBase\_setup.m
    - '1' → plot all except the 3<sup>rd</sup> band defined in DR\_PBase\_setup.m
    - nor: how to normalize spectra. Currently 1
    - plot\_clust: the cluster (defined in DR\_PBase\_setup.m) to be plotted (in matrix form; sequence matter)
  - ◆ Directories: FilePath/
- Outputs
  - ◆ Figure:
    - ICs(plot\_clust)\_(powinc)\_(freq\_name)(nor)(rj)(SD)
    - .fig. Print the png file in PC
    - Note: plot\_clust is hexadecimal

